

Osteoblast Differentiation and Mineralization

Application Note

Background

Osteoblasts (HOB) are specialized fibroblasts that secrete and mineralize the bone matrix. They develop from mesenchymal precursors. The mineralized extracellular matrix is mainly composed of type I collagen and smaller but significant amounts of osteocalcin (OC), matrix gla protein, osteopontin (OPN), bone sialoprotein (BSP), BMPs, TGF- β , and the inorganic mineral hydroxylapatite.

Osteoblast differentiation *in vitro* and *in vivo* can be characterized in three stages: (a) cell proliferation, (b) matrix maturation, and (c) matrix mineralization [1].

In vitro, matrix maturation and mineralization are usually enhanced by growing the cells to complete confluency and by adding specific osteogenic factors [2].


(a) During proliferation, several extracellular matrix proteins (procollagen I, TGF- β , and fibronectin) can be detected. The matrix maturation phase (b) is characterized by maximal expression of alkaline phosphatase (AP).

Finally, at the beginning of matrix

mineralization (c), genes for proteins such as OC, BSP, and OPN are expressed and once mineralization is completed, calcium deposition can be visualized using adequate staining methods.

Analysis of bone cell-specific markers like AP, OC, and collagen type I or detection of functional mineralization is frequently used to characterize osteoblasts *in vitro* [2].

The mineralization process of osteoblasts in *in vitro* culture has also been used as a model for testing the effects of drug treatments and mechanical loading on bone cell differentiation and bone formation [3, 4].



Important: Do not let the cells dry for longer than 30 sec. throughout the entire staining procedure!

Detection of Alkaline Phosphatase*

Proliferating Osteoblasts show alkaline phosphatase (AP) activity, which is greatly enhanced during *in vitro* bone formation. AP activity is therefore a feasible marker for HOB. AP can easily be detected using BCIP/NBT (5-Bromo-4-chloro-3-indolyl phosphate/Nitro blue tetrazolium) as a substrate, which stains cells blue-violet when AP is present.

1. Prepare solutions and reagents

Obtain Saccomanno Fixation Solution (Morphisto, #13881.00250). Dissolve one BCIP/NBT tablet (SigmaFast™ BCIP-NBT; Sigma Aldrich) in 10 ml distilled water to prepare the substrate solution. Store in the dark and use within 2 hours.

Add 0.05% Tween 20 to PBS, w/o Ca⁺⁺/ Mg⁺⁺ (Cat. No. C-40232) to prepare the washing buffer.

2. Wash the cells

Remove the cells from the incubator and carefully aspirate the medium. Carefully wash the cells with PBS.

Note: Do not disrupt the cell monolayer!

3. Fixation of the cells

Carefully aspirate the PBS and add enough Saccomanno Fixation Solution to cover the cellular monolayer. After 60 - 90 seconds gently aspirate the fixation solution and wash the cells with washing buffer.

Note: Longer fixation will lead to irreversible inactivation of AP.

4. Stain the cells

Carefully aspirate the washing buffer and add enough BCIP/NBT substrate solution to cover the cellular monolayer. Incubate at room temperature in the dark for 5-10 min. Check staining progress every 2-3 min.

5. Wash the cells

Carefully aspirate the substrate solution and wash the cell monolayer with washing buffer. Carefully aspirate the washing buffer and add PBS.

6. Analyze the cells

Evaluate staining results. Refer to Fig. 1 for an example of AP detection.




Fig. 1: Alkaline Phosphatase detection

HUVEC (AP negative, upper row) are colorless or faintly bluish, whereas osteoblasts (AP positive, lower row) are dark blue-violet. The higher the AP activity, the more intense the color.

Detection of Alkaline Phosphatase

Please follow the recommended safety precautions for the chemicals used in this procedure!

* AP activity is not limited to osteoblasts. Therefore a second confirmation, e.g. direct staining of extracellular calcium deposits (mineralization), may be necessary.

Use aseptic techniques and a laminar flow bench.

Osteoblast Mineralization

1. Seed Osteoblasts (HOB) on collagen precoated culture plates

Plate 3×10^4 HOB per well on a collagen I coated 24-well tissue culture plate (Corning, Cat. No. #354408). Work in duplicate. Use HOB Growth Medium (C-27001) for one well as a negative control and Osteoblast Mineralization Medium (C-27020) for the other well.

2. Differentiation culture of induced Osteoblasts

Incubate the cells for 17–21 days. Change the medium every third day. Be careful not to disturb the cell monolayer.

Proceed with the protocol in the following section “Detection of Calcium Deposits (Mineralization).”




Fig. 2: Microscopic appearance of HOB after mineralization *in vitro*. Left: Mineralized osteoblasts in Osteoblast Mineralization Medium show vast extracellular calcium deposits, stained in bright orange-red. Right: The negative control in HOB Growth Medium is slightly reddish.




Fig. 3: Macroscopic appearance of HOB after mineralization *in vitro*. The negative control in Osteoblast Growth Medium (upper row) is slightly reddish, whereas the mineralized osteoblasts in Osteoblast Mineralization Medium show vast extracellular calcium deposits, stained in bright orange-red (lower row).

Osteoblast Mineralization

Important: Do not let the cells dry for longer than 30 sec. throughout the entire staining procedure!

Detection of Calcium Deposits (Mineralization)

Osteoblasts can be induced to produce vast extracellular calcium deposits *in vitro*. This process is called mineralization. Calcium deposits are an indication of successful *in vitro* bone formation and can specifically be stained bright orange-red using Alizarin Red S.

1. Prepare solutions and buffers

Use Saccomanno Fixation Solution (Morphisto, #13881.00250). To prepare the Alizarin Red S staining solution dissolve 2 g Alizarin Red S in 90 ml distilled water, mix and adjust the pH to 4.1–4.3 with hydrochloric acid, as necessary. Then, bring up to a final volume of 100 ml with distilled water and filter the dark-brown solution. Store in the dark at 2–8 °C.

Note: The correct pH of the solution is critical. Check pH (at ambient temperature) if the solution is more than 1 month old.

2. Wash the cells

Remove the cells from the incubator and carefully aspirate the medium. Carefully wash the cells with Dulbecco's PBS, w/o Ca⁺⁺/ Mg⁺⁺ (C-40232).

Note: Do not disrupt the cell monolayer!

3. Fixation of the cells

Carefully aspirate the PBS and add enough Saccomanno Fixation Solution to cover the cellular monolayer. After at least 60 min gently aspirate the fixation solution and wash the cells with distilled water.

4. Stain the cells

Immediately before use, pass the required amount of Alizarin Red S staining solution through a 0.22 µm syringe filter equipped with a PES-membrane.

Carefully aspirate the distilled water and add enough filtered Alizarin Red S staining solution to cover the cellular monolayer. Incubate at room temperature in the dark for 45 min.

5. Wash the cells

Carefully aspirate the Alizarin Red S staining solution and wash the cell monolayer four times with 1 ml distilled water. Carefully aspirate the distilled water and add PBS.

6. Analyze the cells

Analyze the sample immediately, as the dye may bleed upon prolonged storage without embedding. Undifferentiated HOB (without extracellular calcium deposits) are slightly reddish, whereas mineralized osteoblasts (with extracellular calcium deposits) are bright orange-red. Refer to Fig. 2 and 3 for an example of osteoblast mineralization.

Detection of Calcium Deposits

Please follow the recommended safety precautions for the chemicals used in this procedure!



References

- [1] Stein GS and Lian JB. Molecular mechanisms mediating developmental and hormone-regulated expression of genes in osteoblasts: an integrated relationship of cell growth and differentiation. In: Noda M, editor. Cellular and molecular biology of bone. Tokyo: Academic Press. p 47–95, 1993.
- [2] Kasperk C. *et al.* Human bone cell phenotypes differ depending on their skeletal site of origin J Clin Endocrinol Metab. Aug;80(8): 2511-7, 1995
- [3] Kostenuik, P.J. *et al.* Skeletal unloading inhibits the in vitro proliferation and differentiation of rat osteoprogenitor cells. Am. J. Physiol. 273, E1133, 1997.
- [4] Kostenuik, P.J. *et al.* Skeletal unloading causes resistance of osteoprogenitor cells to parathyroid hormone and to insulin-like growth factor-I. J. Bone Miner. Res. 14, 21, 1999.

Related Products

Product	Size	Catalog Number
Human Osteoblasts (HOB)	500,000 cryopreserved cells	C-12720
	500,000 proliferating cells	C-12760
Osteoblast Growth Medium (Ready-to-use)	500 ml	C-27001
Osteoblast Mineralization Medium (Ready-to-use)	100 ml	C-27020
DetachKit	30 ml	C-41200
	125 ml	C-41210
	250 ml	C-41220
Cryo-SFM	30 ml	C-29910
	125 ml	C-29912
Dulbecco's PBS, w/o Ca ⁺⁺ / Mg ⁺⁺	500 ml	C-40232
HOB Pellet	> 1 million cells per pellet	C-14071

PromoCell GmbH

Sickingenstr. 63/65
69126 Heidelberg
Germany

Email: info@promocell.com
www.promocell.com

USA/Canada

Phone: 1 – 866 – 251 – 2860 (toll free)
Fax: 1 – 866 – 827 – 9219 (toll free)

Deutschland

Telefon: 0800 – 776 66 23 (gebührenfrei)
Fax: 0800 – 100 83 06 (gebührenfrei)

France

Téléphone: 0800 90 93 32 (ligne verte)
Téléfax: 0800 90 27 36 (ligne verte)

United Kingdom

Phone: 0800 – 96 03 33 (toll free)
Fax: 0800 – 169 85 54 (toll free)

Other Countries

Phone: +49 6221 – 649 34 0
Fax: +49 6221 – 649 34 40

© PromoCell GmbH