

## Specifications

| Product  | Recommended Culture Media* | Plating Density                           | Passage after Thawing | Marker                   | Population Doublings |
|--|----------------------------|---|-----------------------|--------------------------|----------------------|
| Human Nasal Epithelial Cells (HNEpC)           | C-21060                    | 10,000 - 15,000 cells per cm <sup>2</sup> | P2                    | Cytokeratin <sup>+</sup> | > 15                 |
| Human Tracheal Epithelial Cells (HTEpC)        | C-21060                    | 10,000 - 15,000 cells per cm <sup>2</sup> | P2                    | Cytokeratin <sup>+</sup> | > 15                 |
| Human Bronchial Epithelial Cells (HBEpC)       | C-21060                    | 10,000 - 15,000 cells per cm <sup>2</sup> | P2                    | Cytokeratin <sup>+</sup> | > 15                 |
| Human Small Airway Epithelial Cells (HSAEpC)   | C-21070                    | 10,000 - 15,000 cells per cm <sup>2</sup> | P2                    | Cytokeratin <sup>+</sup> | > 15                 |
| Human Renal Epithelial Cells (HREpC)           | C-26001<br>C-26030         | 10,000 - 15,000 cells per cm <sup>2</sup> | P2                    | Cytokeratin <sup>+</sup> | > 15                 |
| Human Renal Cortical Epithelial Cells (HRCepC) | C-26001<br>C-26030         | 10,000 - 15,000 cells per cm <sup>2</sup> | P2                    | Cytokeratin <sup>+</sup> | > 15                 |
| Human Placental Epithelial Cells (HPIEpC)      | C-26040                    | 10,000 - 15,000 cells per cm <sup>2</sup> | P2                    | Cytokeratin <sup>+</sup> | > 5                  |

## Related Products

| Product  | Size       | Catalog Number |
|--|------------|----------------|
| Airway Epithelial Cell Growth Medium (Ready-to-use)        | 500 ml     | C-21060        |
| Airway Epithelial Cell Growth Medium Kit                   | 500 ml     | C-21160        |
| Airway Epithelial Cell Basal Medium                        | 500 ml     | C-21260        |
| Airway Epithelial Cell Basal Medium, phenol red-free       | 500 ml     | C-21265        |
| Airway Epithelial Cell Growth Medium SupplementMix         | for 500 ml | C-39165        |
| Airway Epithelial Cell Growth Medium SupplementPack        | for 500 ml | C-39160        |
| Small Airway Epithelial Cell Growth Medium (Ready-to-use)  | 500 ml     | C-21070        |
| Small Airway Epithelial Cell Growth Medium Kit             | 500 ml     | C-21170        |
| Small Airway Epithelial Cell Basal Medium                  | 500 ml     | C-21270        |
| Small Airway Epithelial Cell Basal Medium, phenol red-free | 500 ml     | C-21275        |
| Small Airway Epithelial Cell Growth Medium SupplementMix   | for 500 ml | C-39175        |
| Small Airway Epithelial Cell Growth Medium SupplementPack  | for 500 ml | C-39170        |
| Mammary Epithelial Cell Growth Medium (Ready-to-use)       | 500 ml     | C-21010        |
| Mammary Epithelial Cell Growth Medium Kit                  | 500 ml     | C-21110        |
| Mammary Epithelial Cell Basal Medium                       | 500 ml     | C-21210        |
| Mammary Epithelial Cell Basal Medium, phenol red-free      | 500 ml     | C-21215        |
| Mammary Epithelial Cell Growth Medium SupplementMix        | for 500 ml | C-39115        |
| Mammary Epithelial Cell Growth Medium SupplementPack       | for 500 ml | C-39110        |

\*The catalog numbers in this table are for media in ready-to-use packaging.

| Product   | Size                       | Catalog Number                |
|---|----------------------------|-------------------------------|
| Renal Epithelial Cell Growth Medium (Ready-to-use)      | 500 ml                     | C-26001                       |
| Renal Epithelial Cell Growth Medium Kit                 | 500 ml                     | C-26011                       |
| Renal Epithelial Cell Basal Medium                      | 500 ml                     | C-26021                       |
| Renal Epithelial Cell Basal Medium, phenol red-free     | 500 ml                     | C-26022                       |
| Renal Epithelial Cell Growth Medium SupplementMix       | for 500 ml                 | C-39604                       |
| Renal Epithelial Cell Growth Medium SupplementPack      | for 500 ml                 | C-39600                       |
| Renal Epithelial Cell Growth Medium 2 (Ready-to-use)    | 500 ml                     | C-26030                       |
| Renal Epithelial Cell Growth Medium 2 Kit               | 500 ml                     | C-26130                       |
| Renal Epithelial Cell Basal Medium 2                    | 500 ml                     | C-26230                       |
| Renal Epithelial Cell Basal Medium 2, phenol red-free   | 500 ml                     | C-26235                       |
| Renal Epithelial Cell Growth Medium 2 SupplementMix     | for 500 ml                 | C-39606                       |
| Renal Epithelial Cell Growth Medium 2 SupplementPack    | for 500 ml                 | C-39605                       |
| Placental Epithelial Cell Growth Medium (Ready-to-use)  | 500 ml                     | C-26040                       |
| Placental Epithelial Cell Growth Medium Kit             | 500 ml                     | C-26140                       |
| Placental Epithelial Cell Basal Medium                  | 500 ml                     | C-26240                       |
| Placental Epithelial Cell Basal Medium, phenol red-free | 500 ml                     | C-26245                       |
| Placental Epithelial Cell Growth Medium SupplementMix   | for 500 ml                 | C-39612                       |
| Placental Epithelial Cell Growth Medium SupplementPack  | for 500 ml                 | C-39610                       |
| DetachKit   | 30 ml<br>125 ml<br>250 ml  | C-41200<br>C-41210<br>C-41220 |
| Cryo-SFM  | 30 ml<br>125 ml            | C-29910<br>C-29912            |
| HNEpC Pellet  | 1 million cells per pellet | C-14062                       |
| HTepC Pellet  | 1 million cells per pellet | C-14064                       |
| HBEpC Pellet  | 1 million cells per pellet | C-14063                       |
| HSAEpC Pellet   | 1 million cells per pellet | C-14065                       |
| HREpC Pellet  | 1 million cells per pellet | C-14067                       |
| HRCEpC Pellet   | 1 million cells per pellet | C-14068                       |
| HPLEpC Pellet   | 1 million cells per pellet | C-14069                       |
| BPE-15  | 2 ml                       | C-30020                       |
| BPE-26  | 2 ml                       | C-30021                       |
| EGF, human, recombinant                                 | 500 µg                     | C-60170                       |

**PromoCell GmbH**

Sickingenstr. 63/65  
69126 Heidelberg  
Germany

Email: [info@promocell.com](mailto:info@promocell.com)  
[www.promocell.com](http://www.promocell.com)

**USA/Canada**

Phone: 1 – 866 – 251 – 2860 (toll free)  
Fax: 1 – 866 – 827 – 9219 (toll free)

**Deutschland**

Telefon: 0800 – 776 66 23 (gebührenfrei)  
Fax: 0800 – 100 83 06 (gebührenfrei)

**France**

Téléphone: 0800 90 93 32 (ligne verte)  
Téléfax: 0800 90 27 36 (ligne verte)

**United Kingdom**

Phone: 0800 – 96 03 33 (toll free)  
Fax: 0800 – 169 85 54 (toll free)

**Other Countries**

Phone: +49 6221 – 649 34 0  
Fax: +49 6221 – 649 34 40

© PromoCell GmbH

## Instruction Manual

| Product  | Size                        | Catalog Number |
|--|-----------------------------|----------------|
| Human Nasal Epithelial Cells (HNEpC)           | 500,000 cryopreserved cells | C-12620        |
|  | 500,000 proliferating cells | C-12621        |
| Human Tracheal Epithelial Cells (HTEpC)        | 500,000 cryopreserved cells | C-12644        |
|  | 500,000 proliferating cells | C-12645        |
| Human Bronchial Epithelial Cells (HBEpC)       | 500,000 cryopreserved cells | C-12640        |
|  | 500,000 proliferating cells | C-12641        |
| Human Small Airway Epithelial Cells (HSAEpC)   | 500,000 cryopreserved cells | C-12642        |
|  | 500,000 proliferating cells | C-12643        |
| Human Renal Epithelial Cells (HREpC)           | 500,000 cryopreserved cells | C-12665        |
|  | 500,000 proliferating cells | C-12666        |
| Human Renal Cortical Epithelial Cells (HRCEpC) | 500,000 cryopreserved cells | C-12660        |
|  | 500,000 proliferating cells | C-12662        |
| Human Placental Epithelial Cells (HPIEpC)      | 500,000 cryopreserved cells | C-12395        |
|  | 500,000 proliferating cells | C-12397        |

### Product Description

The epithelium is the interface between the body and the external environment and covers all exterior surfaces and interior lumina. Depending on the tissue of origin, the functions of epithelial cells are diverse and include absorption, secretion, protection, transcellular transport, and sensation.

PromoCell offers a range of Epithelial Cells produced at PromoCell's cell culture facility. The cells are isolated from normal human adult tissues of the airways (nasal mucosa, trachea, bronchi, and distal respiratory tract), kidney (whole kidney and renal cortex), and the placenta (amniotic membrane). Shortly after isolation, all PromoCell Human Epithelial Cells are cryopreserved at passage 2 (P2) using PromoCell's proprietary, serum-free freezing medium, Cryo-SFM. Each cryo vial contains more than 500,000 viable cells after thawing. Proliferating cell cultures are made

from cryopreserved cells that have been thawed and cultured for three days at PromoCell.

### Quality Control

Rigid quality control tests are performed for each lot of Epithelial Cells.

They are tested for cell morphology, adherence rate and cell viability. Furthermore, flow cytometric analyses for the cell-type specific marker cytokeratin are carried out for each lot (see page 5). Growth performance is tested through multiple passages under culture conditions without antibiotics and antimycotics. In addition, all cells have been tested for the absence of HIV-1, HIV-2, HBV, HCV, HTLV-1, HTLV-2 and microbial contaminants (fungi, bacteria, and mycoplasma). A detailed certificate of analysis (CoA) for each lot can be downloaded at: [www.promocell.com/coa](http://www.promocell.com/coa)

### Intended Use

PromoCell Epithelial Cells are for *in vitro* research use only and not for diagnostic or therapeutic procedures.

### Warning

Although tested negative for HIV-1, HIV-2, HBV, HCV, HTLV-1 and HTLV-2, the cells – like all products of human origin – should be handled as potentially infectious. No test procedure can completely guarantee the absence of infectious agents.

**Follow appropriate safety precautions!**

*After delivery, cryopreserved cells should be stored in liquid nitrogen or seeded directly (see page 2). Proliferating cells have to be processed immediately (see page 3).*

*Use aseptic techniques and a laminar flow bench.*

### Protocol for Cryopreserved Cells

**Straight after arrival, store the cryopreserved cells in liquid nitrogen or seed them immediately.**

**Note:** Storage at  $-80^{\circ}\text{C}$  is not sufficient for cell preservation and causes irreversible cell damage.

#### 1. Prepare the medium

Calculate the needed culture surface area according to the plating density (see page 5) and the lot-specific cell numbers stated on the certificate of analysis. Fill the appropriate volume of PromoCell Growth Medium (at least 9 ml per vial of cells) in cell culture vessels. Place the vessels in an incubator ( $37^{\circ}\text{C}$ , 5%  $\text{CO}_2$ ) for 30 minutes.



#### 2. Thaw the cells

Remove the cryovial from the liquid nitrogen container and immediately place it on dry ice – even for short transportation. Under a laminar flow bench, briefly twist the cap a quarter turn to relieve pressure, then re-tighten. Immerse the vial into a water bath ( $37^{\circ}\text{C}$ ) just up to the screw cap for 2 minutes. Ensure that no water enters the thread of the screw cap.



#### 3. Disinfect the vial and seed the cells

Thoroughly rinse the cryovial with 70% ethanol under a laminar flow bench. Then, aspirate the excess ethanol from the thread area of the screw cap. Open the vial and transfer the cells to a cell culture vessel containing the prewarmed medium from step 1.



#### 4. Incubate the cells

Place the vessel in an incubator ( $37^{\circ}\text{C}$ , 5%  $\text{CO}_2$ ) for cell attachment. Replace the medium after 16–24 hours and every two to three days thereafter. The cells should be subcultured, according to the subcultivation protocol (see page 4), once they have reached 70–90% confluency.



*Start immediately after delivery.  
Use aseptic techniques and a laminar flow bench.*

### Protocol for Proliferating Cells

#### 1. Incubate the cells

Unpack the culture vessel, do not open the cap, and immediately place it in an incubator (37°C, 5% CO<sub>2</sub>) for 3 hours to allow the cells to recover from the transportation.




#### 2. Replace the transport medium

Carefully open the vessel, rinse the inner side of the cap with 70% ethanol, and let air dry. Aspirate the transport medium from the vessel. Add 10 ml of the appropriate PromoCell Cell Growth Medium.



#### 3. Check and incubate the cells

Check the cell density. Open the lid half a turn and place the vessel in an incubator (37°C, 5% CO<sub>2</sub>). Change the medium every two to three days. The cells should be subcultured, according to the subcultivation protocol (see page 4), once they have reached > 70% confluency.



*Use aseptic techniques and a laminar flow bench.*

## Subcultivation Protocol

### 1. Prepare the reagents and wash the cells

Place the PromoCell DetachKit at room temperature for at least 30 minutes to adjust the temperature of the reagents. Carefully aspirate the medium from the culture vessel. Add 100  $\mu$ l Hepes BSS Solution per  $\text{cm}^2$  of vessel surface to wash the cells and agitate the vessel carefully for 15 seconds.



### 2. Detach the cells

Carefully aspirate the Hepes BSS from the culture vessel. Add 100  $\mu$ l Trypsin/EDTA Solution per  $\text{cm}^2$  of vessel surface. **Note:** We recommend detaching the cells at room temperature. Close the vessel and examine the cells under a microscope. When the cells start to detach, gently tap the side of the vessel to loosen the remaining cells.



### 3. Neutralize the trypsin and harvest the cells

Add 100  $\mu$ l Trypsin Neutralization Solution per  $\text{cm}^2$  of vessel surface and gently agitate. Carefully aspirate the cell suspension and transfer it to a centrifugation tube. Spin down the cells for 3 minutes at 220 x g.



### 4. Incubate the cells

Discard the supernatant (step 1), add 1 ml of the appropriate PromoCell Cell Growth Medium (step 2), and resuspend the cells by carefully pipetting up and down. Plate the cells according to the recommended seeding density in new cell culture vessels containing prewarmed PromoCell Cell Growth Medium. Place the vessels in an incubator (37°C, 5%  $\text{CO}_2$ ) and change the media every two to three days.

