

Standardized culture of cancer cell lines under serum- and xeno-free conditions

Application note

This application note describes the switch of an existing cancer cell culture to our standardized Cancer Cell Line Medium XF and subsequent routine culture.

Traditionally, established cell lines have been propagated almost exclusively in standard culture media supplemented with significant amounts (2–20%) of fetal calf serum (FCS). Like other undefined media supplements, fetal calf serum has unwanted physiological, genetic and epigenetic cellular effects [1–4] and is known to cause enormous experimental variability and distort readouts and results, e.g., in drug screenings and

hormone-related studies [5]. With the aim of achieving a universal, consistent and xeno-free environment for culturing most of the commonly used human cancer cell lines, our Cancer Cell Line Medium XF was designed as a defined formulation with the exception of highly purified human plasma-derived albumin. It has no ill-defined components such as fetal calf serum, extracts or hydrolysates and exhibits very low lot-to-lot variability. Due to the flexible nature of our Cancer Cell Line Medium XF, you can use it to support a broad range of cancer cell cells from tumor biopsies while maintaining the tumor-initiating cancer stem cells (CSCs).

It further allows you to cultivate or co-cultivate a range of different cancer cells, cancer-associated cells, and non-cancer cells such as immune or stromal cells.

Since the formulation does not contain any cellular attachment factors, it is necessary to coat culture vessels with fibronectin, vitronectin or other suitable attachment substrates. Being broadly usable across all common adherently growing cancer cell lines, this new culture medium is a cost-effective solution for ensuring efficient, genuinely standardized routine cultures.




Fig. 1A: Expansion performance of the HT1080 fibrosarcoma cell line in the Cancer Cell Line Medium XF compared with conventional culture conditions. HT1080 cells were plated with 5,000 cells/cm² in Cancer Cell Line Medium XF on fibronectin-coated vessels (red) or in DMEM+ 2 mM L-Glutamine + 10% FCS (grey). Subsequently, the cells were cultured for 10 consecutive passages with a passage interval of 3 to 4 days.

B,C: Morphology of HT1080 fibrosarcoma cells cultured in Cancer Cell Line Medium XF. Exemplary images of HT1080 on day three after subculture (P7) are shown in the Cancer Cell Line Medium XF (B) and conventional culture conditions (C) (100x magnification).

Background

For cost reasons, immortal cancer cell lines are widely used in research instead of primary cells. They are accepted, well-defined model systems that serve as a constant source of cells while avoiding the limitations posed by the finite lifespan of normal primary cells [6–9]. As long as the limitations inherent in substituting cell lines for primary cells are considered, tumor cells can be used in *in vitro* models to reflect certain aspects and functionalities of terminally differentiated primary human cells [7]. However, poorly de-

fined culture media components such as fetal calf serum are a well-known and significant source of variations [1–5]. These can cancel out the advantages of using cancer cell lines as a uniform source of different types of cells while endangering the reliability of experimental results obtained with them. As it has been demonstrated for hormone responsive cell lines, immunosystem-related cell types and stem cells, undefined conventional culture media can interfere with cell properties and experimental readouts, for instance by falsi-

fying their responses to drugs, provoking un-specific bogus immune responses or causing unwanted differentiation of stem cells [1, 3, 10]. Cells can also be significantly altered by culture conditions of this kind [11]. It follows that a controlled culture environment is key for obtaining more accurate results with cell line models and facilitates data analysis and interpretation [4], thus taking a major step forward in exploiting the unique features of permanent cell lines as cost-effective *in vitro* research models.




Fig. 2: Morphology of MCF-7 breast carcinoma cells in the Cancer Cell Line Medium XF. MCF-7 cells were plated at 10,000 cells per cm^2 in fibronectin-coated vessels and cultured for three passages in the Cancer Cell Line Medium XF. The cells exhibited efficient proliferation as well as a typical – but slightly more compact – morphology as compared to traditional culture media (not shown). The image was taken 5 days after seeding at 100x magnification.

Routine cancer cell culture protocol

This protocol describes how a human cancer cell line can be switched from a standard medium to the Cancer Cell Line Medium XF.

I. Switching an existing cell culture to the Cancer Cell Line Medium XF

Materials

- Proliferating culture of a human cancer cell line
- Cancer Cell Line Medium XF
- Adhesion factors: Human fibronectin (C-43060) or vitronectin
- Phosphate buffered saline (PBS) without Ca^{2+} / Mg^{2+}
- Accutase (C-41310) or DetachKit (C-41210)
- Tissue-culture-treated cell culture vessels

Use aseptic techniques and a laminar flow bench.

1

Coat the culture vessel

The serum- and xeno-free formulation does not contain attachment factors. Thus coating of the surface of the cell culture vessel with an appropriate adhesion factor is usually needed. Table 1 (p. 5) shows an overview of cell lines and surface coatings tested with the Cancer Cell Line Medium XF.

For the establishment of the culture conditions, it is recommended to test fibronectin and vitronectin coating: Coat the culture vessel

with 10 $\mu\text{g}/\text{ml}$ human (or bovine) fibronectin or 5 $\mu\text{g}/\text{ml}$ vitronectin according to the instruction manual of the product. Use 100 μl of diluted coating solution per cm^2 of culture surface. (Final concentration: fibronectin 1 $\mu\text{g}/\text{cm}^2$ and vitronectin 0.5 $\mu\text{g}/\text{cm}^2$).

Note: If not used immediately, the sealed vessel may be stored for up to 3 months at 2–8°C for later use.

2

Harvest cells from your existing culture

Harvest and count cells from an established culture of the appropriate cell line using your standard method. Resuspend them in Cancer Cell Line Medium XF.

3

Plate the cells

Plate the cells at a density of 5,000–10,000 cells/ cm^2 . When seeding the cells for the first time in the Cancer Cell Line Medium XF, use approximately 200 μl of medium per cm^2 of culture surface, e.g., 5 ml for a T25 flask.

4

Let the cells grow

Incubate the plated cells at 37°C and 5% CO_2 . Change the medium every 2 to 3 days.

Note: Adaption of cell cultures to the Cancer Cell Line Medium XF is not required. With some cell lines, proliferation may be somewhat reduced after initiating the culture but this should normalize after one to three passages.

Subculture the cells

Once the cells have reached 70–80% confluence, wash the culture twice with ambient tempered PBS without $\text{Ca}^{2+}/\text{Mg}^{2+}$ and then incubate the cells for 5–10 minutes with $150 \mu\text{l}/\text{cm}^2$ Accutase at 37°C . After the first 5 minutes of incubation, monitor the detachment process visually. When the cells start to detach, facilitate their complete dislodgement by tapping the flask. Add the same volume of Cancer Cell Line Medium XF to the detached cells and spin down for 5 minutes at $300 \times g$ at room temperature. Carefully aspirate the supernatant and gently resuspend the cell pellet in an adequate amount of Cancer Cell Line Medium XF. Seed the cells into new fibronectin-coated vessels and incubate them further at 37°C and 5% CO_2 . Use approx. $300\text{--}400 \mu\text{l}$ of medium per cm^2 of culture surface for the subsequent cultivation. Continue incubation of the cultures at 37°C and 5% CO_2 .

II. Routine culture using the Cancer Cell Line Medium XF

Materials

- Proliferating culture of a human cancer cell line
- Cancer Cell Line Medium XF
- Adhesion factors: Human fibronectin (C-43060) or vitronectin
- Phosphate buffered saline (PBS) without $\text{Ca}^{2+}/\text{Mg}^{2+}$ (C-40232)
- Accutase (C-41310) or DetachKit (C-41210)
- Tissue-culture-treated cell culture vessels

Use aseptic techniques and a laminar flow bench.

This protocol describes the routine culture of a human cancer cell line already transferred to the Cancer Cell Line Medium XF.

1

Coat the culture vessel

Coat the culture vessel with the suitable adhesion factors as previously tested (see page 3).

2a

Harvest cells from your existing culture

Harvest and count cells from an established culture of the appropriate cell line using your standard method. Resuspend them in Cancer Cell Line Medium XF.

2b

Thaw cells from cryo-stock

Cryopreserved cells cultured previously in Cancer Cell Line Medium XF may also be thawed directly using this medium.

3

Plate the cells

Plate the cells at the appropriate density ($5,000\text{--}10,000 \text{ cells}/\text{cm}^2$). Use approximately $300\text{--}400 \mu\text{l}$ of Cancer Cell Line Medium XF per cm^2 of culture surface, e.g., 7.5–10 ml for a T25 flask.

4

Let the cells grow

Incubate the plated cells at 37°C and 5% CO_2 . Change the medium every 2 to 3 days.

Subculture the cells

Once the cells have reached 70–80% confluence, wash the culture twice with ambient tempered PBS without Ca²⁺/Mg²⁺ and then incubate the cells for 5–10 minutes with 150 µl/cm² Accutase at 37°C. After the first 5 minutes of incubation, monitor the detachment process visually. When the cells start to detach, facilitate their complete dislodgement by tapping the flask. Add the same volume of Cancer Cell Line Medium XF to the detached cells and spin down for 5 minutes at 300 x g at room temperature. Carefully aspirate the supernatant and gently resuspend the cell pellet in an adequate amount of Cancer Cell Line Medium XF. Seed the cells into new fibronectin-coated vessels and incubate them further at 37°C and 5% CO₂. Use approx. 300–400 µl of medium per cm² of culture surface for the subsequent cultivation. Continue incubation of the cultures at 37°C and 5% CO₂.

Tissue	Tested cell line	Cell line origin	Remarks
Abelson murine leukemia virus-induced tumor	RAW264.7	Mouse Macrophage	Coat with fibronectin: 1 µg/cm ²
Bone marrow	KG-1	Human acute myelogenous leukemia (suspension)	No coating required
Brain	BV2	Immortalized murine primary microglial cells	Coat with fibronectin: 1 µg/cm ²
Brain	A172	Human Glioblastoma	Coat with fibronectin: 1 µg/cm ²
Brain, neuroectodermal	Neuro-2a	Mouse Neuroblastoma	Coat with fibronectin: 1 µg/cm ²
Breast	MCF-7	Pleural effusion of metastatic human breast adenocarcinoma	Coat with fibronectin: 1 µg/cm ²
Breast	HCC38	Human breast ductal carcinoma cell line	Coat with collagen 6 µg/cm ²
Breast	HCC1143	Human ductal carcinoma cell line	Coat with collagen 6 µg/cm ²
Breast	HCC1937	Human ductal carcinoma cell line	Coat with collagen 6 µg/cm ²
Colon	HT-29	Human colon adenocarcinoma	Coat with vitronectin: 0.5 µg/cm ²
Connective tissue	HT 1080	Human fibrosarcoma	Coat with fibronectin: 1 µg/cm ²
Liver	HepG2	Hepatocellular carcinoma of the human liver	Coat with vitronectin: 0.5 µg/cm ²
Lung	A-549	Human lung carcinoma	Coat with vitronectin: 0.5 µg/cm ²
Peripheral blood	BDCM	B lymphoblast cell line	No coating required
Peripheral blood	HCC38 BL	EBV transformed B lymphoblast cell line	No coating required
Peripheral blood	MV-4-11	Human acute myelogenous leukemia (suspension)	No coating required
Prostate	LNCaP	Lymph node metastasis of human prostate adenocarcinoma	3D culture in 3D Tumorsphere Growth Medium XF (C-28070) is recommended
Skeletal muscle	C2C12	Mouse Myoblasts	Coat with fibronectin: 1 µg/cm ²
Skin	B16-F10	Mouse Melanoma	Coat with fibronectin: 1 µg/cm ²

Table 1: Overview of cell lines and surface coatings tested with the Cancer Cell Line Medium X. Cells were seeded at a density of 10,000 cells/cm².

Products

Media	Size	Catalog number
Cancer Cell Line Medium XF	250 ml	C-28077
Accutase-Solution	100 ml	C-41310

Related Products

Media	Size	Catalog number
3D Tumorsphere Medium XF	250 ml	C-28070
Primary Cancer Culture System	250 ml	C-28081
Dulbecco's PBS, without Ca ⁺⁺ /Mg ⁺⁺	500 ml	C-40232
DetachKit	3 × 125 ml	C-41210
Cryo-SFM	30 ml / 125 ml	C-29910/C-29912

References

1. Sulit, H.L., et al., Human tumor cells grown in fetal calf serum and human serum: influences on the tests for lymphocyte cytotoxicity, serum blocking and serum arming effects. *Int J Cancer*, 1976. 17(4): pp. 461–8.
2. Darro, F., et al., Characterization of the differentiation of human colorectal cancer cell lines by means of Voronoi diagrams. *Cytometry*, 1993. 14(7): pp. 783–92.
3. Sedelaar, J.P.M. and J.T. Isaacs, Tissue Culture Media Supplemented with 10% Fetal Calf Serum Contains a Castrate level of Testosterone. *The Prostate*, 2009. 69(16): pp. 1724–1729.
4. Baker, M., Reproducibility: Respect your cells! *Nature*, 2016. 537(7620): pp. 433–435.
5. Sikora, M.J., et al., Endocrine Response Phenotypes Are Altered by Charcoal-Stripped Serum Variability. *Endocrinology*, 2016. 157(10): pp. 3760–3766.
6. Kim, H.S., Y.J. Sung, and S. Paik, Cancer Cell Line Panels Empower Genomics-Based Discovery of Precision Cancer Medicine. *Yonsei Med J*, 2015. 56(5): pp. 1186–98.
7. Masters, J.R., Human cancer cell lines: fact and fantasy. *Nat Rev Mol Cell Biol*, 2000. 1(3): pp. 233–6.
8. Kondo, T., Stem cell-like cancer cells in cancer cell lines. *Inflammation and Regeneration*, 2007. 27(5): pp. 506–511.
9. Fillmore, C.M. and C. Kuperwasser, Human breast cancer cell lines contain stem-like cells that self-renew, give rise to phenotypically diverse progeny and survive chemotherapy. *Breast Cancer Res*, 2008. 10(2): pp. R25.
10. Hurt, E.M., et al., Identification of Vitronectin as an Extrinsic Inducer of Cancer Stem Cell Differentiation and Tumor Formation. *Stem Cells (Dayton, Ohio)*, 2010. 28(3): pp. 390–398.
11. Martin, M.J., et al., Human embryonic stem cells express an immunogenic nonhuman sialic acid. *Nat Med*, 2005. 11(2): pp. 228–32.

PromoCell GmbH
Sickingenstr. 63/65
69126 Heidelberg
Germany

USA/Canada
Phone: 1 – 866 – 251 – 2860 (toll free)
Fax: 1 – 866 – 827 – 9219 (toll free)

Deutschland
Telefon: 0800 – 776 66 23 (gebührenfrei)
Fax: 0800 – 100 83 06 (gebührenfrei)

France
Téléphone: 0800 – 90 93 32 (ligne verte)
Téléfax: 0800 – 90 27 36 (ligne verte)

United Kingdom
Phone: 0800 96 03 33 (toll free)
Fax: 0800 169 85 54 (toll free)

Other Countries
Phone: +49 6221 – 649 34 0
Fax: +49 6221 – 649 34 40