

HLA-Typed Human Primary Cells

Discover our unique range of human leukocyte antigen (HLA)-typed cell donors. Choose from our large inventory of HLA-typed donors and start your experiments immediately without wasting time and money on donor screenings.

Musculoskeletal System

Skeletal Muscle Cells
Chondrocytes
Osteoblasts



Cardiovascular System

Smooth Muscle Cells
Endothelial Cells
Cardiac Myocytes
Aortic and Cardiac Fibroblasts



Skin & Connective Tissue

Fibroblasts
Keratinocytes
Melanocytes
Dermal Endothelial Cells
Follicle Dermal Papilla Cells



Respiratory System

Airway Epithelial Cells
Pulmonary Fibroblasts
Pulmonary Endothelial Cells
Smooth Muscle Cells



Kidney

Renal Epithelial Cells



Adipose Tissue

Preadipocytes
Mesenchymal Stem Cells



Bone Marrow

Mesenchymal Stem Cells



Uterus

Smooth Muscle Cells
Fibroblasts
Endothelial Cells




Umbilical Cord

Mesenchymal Stem Cells
Mononuclear Cells
Smooth Muscle Cells
CD34+ Progenitor Cells



Peripheral Blood

Macrophages
Mononuclear Cells



Advantages

- Minimize off-target interactions preclinically by testing translatability in a human cell setting
- Large inventory featuring more than 100 HLA-typed donors in stock and ready to ship
- Complete HLA typing report is delivered for each donor
- 4 digit resolution typing for HLA-A, -B, -C, -DPA1/DPB1, -DQA1/DQB1, and -DRB1/DRB345

Discover

Learn more:

<https://promocell.com/services/hla-typed-human-primary-cells/>

Contact us:

info@promocell.com

PromoCell GmbH

Sickingenstr. 63/65
69126 Heidelberg
Germany

info@promocell.com
www.promocell.com

USA/Canada

Phone: 1 – 866 – 251 – 2860 (toll free)
Fax: 1 – 866 – 827 – 9219 (toll free)

Deutschland

Telefon: 0800 – 776 66 23 (gebührenfrei)
Fax: 0800 – 100 83 06 (gebührenfrei)

France

Téléphone: 0800 – 90 93 32 (ligne verte)
Téléfax: 0800 – 90 27 36 (ligne verte)

United Kingdom

Phone: 0800 96 03 33 (toll free)
Fax: 0800 169 85 54 (toll free)

Other Countries

Phone: +49 6221 – 649 34 0
Fax: +49 6221 – 649 34 40

© PromoCell GmbH